


Animal Care | *Systems*

Handbook for Architects, Engineers, and Building Managers



Contents

- 3 Introduction
- 4 Optimice®
- 5 Optirat® GenII
- 6 Optirat® Plus
- 7 Workspace Requirements
- 8 Airflow Requirements
- 9 Static Pressure Requirements
- 10 Sterilization Requirements
- 11 Exhaust Connections
- 12 Auto-water Systems
- 13 Doorways

Introduction


This guide will assist architecture firms designing buildings tailored to our technology. It provides detailed physical dimensions of the major components of Animal Care Systems racks. It also covers the requirements of air flow, air pressure, and workspace.

Animal Care Systems technology is unique within the animal research industry. Vivaria erected 25 years ago or more were designed around library-style racks with rectangular cages. But our systems are unique.


The slim carousel design fits 30 percent more cages in a given room than other manufacturers' racks, and because Animal Care Systems technology requires no electricity, and animal-generated heat is exhausted immediately, there is no heat load limit for rooms with our racks.



Optimice



The rack frame is 304 stainless steel.

	Optimice Rack					
	Dimensions Inches (cm)				Weight Pounds (kg)	
	Depth	Width	Height	Auto-Water Height	Basic	Auto-Water
Full rack with card holders	40.5 (102.9)	40.5 (102.9)	76.9 (195.3)	77.9 (197.8)	613 (278)	659 (299)
Full rack no card holders	39.8 (101.1)	39.8 (101.1)			604 (274)	650 (295)
Empty rack	33.8 (85.9)	39.8 (101.1)			326 (148)	366 (166)

Racks arrive on standard 48-by-40-inch pallets. Other equipment ships on pallets of up to 52 inches long and 44 inches in wide.

Contact [Animal Care Systems](http://www.animalcaresystems.com) for dimensions of Optimice cage components.



Sample Room Layout

37 feet 9 inches by 30 feet 2 inches


50 Optimice Racks

4 Change Stations


5,000 Mouse Cages

Animal Care Systems provides
Revit files upon request.

Optirat GenII




The rack frame is 304 stainless steel.

	Optirat GenII Rack					
	Dimensions Inches (cm)				Weight Pounds (kg)	
	Depth	Width	Height	Auto-Water Height	Basic	Auto-Water
Full rack	39.4 (100.1)	39.7 (100.8)	76.3 (193.8)	78 (198.1)	534 (242)	524 (238)
Empty rack	33.8 (85.8)	39.7 (100.8)			311 (141)	296 (134)

Racks arrive on standard 48-by-40-inch pallets. Other equipment ships on pallets of up to 52 inches long and 44 inches in wide.


Contact [Animal Care Systems](http://www.animalcaresystems.com) for dimensions of Optirat GenII cage components.



Sample Room Layout
 23 feet 6 inches by 20 feet 2 inches
 20 Optirat GenII Racks
 2 Change Stations
 840 Rat Cages

Animal Care Systems provides
 Revit files upon request.

Optirat Plus




	Optirat Plus Rack					
	Dimensions Inches (cm)				Weight Pounds (kg)	
	Depth	Width	Height	Auto-Water Height	Basic	Auto-Water
Full rack	41.0 (104.1)	39.7 (100.8)			503 (228)	500 (227)
Full rack, carousel rotated for precise entryways	39.4 (1000)	39.7 (100.8)	78 (198.1)	79.5 (201.9)	503 (228)	500 (227)
Empty rack	33.8 (85.8)	39.7 (101)			269 (122)	284 (129)

Racks arrive on standard 48-by-40-inch pallets. Other equipment ships on pallets up to 52 inches long and 44 inches wide.


Contact [Animal Care Systems](http://www.animalcaresystems.com) for dimensions of Optirat Plus cage components.

Animal Care Systems provides Revit files upon request.



Sample Room Layout
20 feet 8 inches by 10 feet 2 inches
10 Optirat Plus Racks
1 Change Station
300 Rat Cages

Workspace Requirements



The rack becomes inaccessible if a stationary object enters a red zone.

Airflow Requirements

Animal Care Systems racks use buildings' HVAC systems to provide negative pressure airflow. Two hoses are recommended, but one hose is viable.


Airflow Requirements for Opti Racks and Cages			
Rack Type	Functional Exhaust Airflow Range, CFM	Anemometer Measurement Range, Per Hose, m/s	
		Two Hoses	One Hose
Optimice	45 - 60	4.0 - 5.25	8.0 - 10.5
Ergo80	38 - 50	3.4 - 4.5	6.8 - 9.0
Ergo70	34 - 45	3.0 - 3.9	6.0 - 7.8
Ergo50	25 - 33	2.2 - 2.9*	4.3 - 5.7
Optirat/GenII/Plus	55 - 82	4.8 - 7.2	11 - 14.4

*A single hose is standard for Ergomice 50

Airflow Requirements for Single-use Caging Within Optimice and Ergomice Racks			
Rack Type	Functional Exhaust Airflow Range, CFM	Anemometer Measurement Range, Per Hose, m/s	
		Two Hoses	One Hose
Optimice	60 - 75	5.3 - 6.6	10.6 - 13.2
Ergo80	51 - 64	4.5 - 5.6	9.0 - 11.2
Ergo70	45 - 56	3.9 - 4.9	7.8 - 9.8
Ergo50	33 - 41	5.7 - 7.1*	11.4 - 14.2

*A single hose is standard for Ergomice 50

Static Pressure Requirements



The dashed box in the line graph represents the 45 to 60 CFM range.

Static pressure was measured in the following locations:

Plenum, 6 in.: Inside the central plenum immediately before the transition to the 6-inch universal box.

Plenum, 4 in.: Inside the central plenum immediately before the transition to the 4-inch universal box.

Hose: Inside the exhaust hoses at the point where they connect to the universal box.

Rack: Inside the exhaust hoses at the point where they connect to the Opti racks.

Sterilization Requirements

Animal Care Systems racks are compatible with steam, dry, and cold sterilization methods.

Contact Animal Care Systems about information on dry heat sterilization options.

<u>Rack Sterilization Temperature Limit</u>
250 F (121 C)

<u>Cage Sterilization Temperature Limit</u>		
	Sterilization Temperature	Maximum Duration
Polycarbonate	up to 250 F (121 C)	20 minutes
Polysulfone	up to 285 F (141 C)	30 minutes



Exhaust Connections

Animal Care Systems racks are compatible with any exhaust connection, including exhaust plates, underceiling structures, and dedicated exhaust drops; however we recommend dedicated exhaust drops.

To achieve the recommended cubic feet per minute, Animal Care Systems recommends ceiling drop diameters of 6 inches or 4 inches; however if two racks are going to be connected to a single drop, we recommend only 6-inch drops instead of 4-inch drops because they impose a smaller load on the HVAC system. The fittings can connect to the building's HVAC system in multiple ways: Universal boxes are used for connecting two racks per drop, while V boxes are recommended for a single rack per drop.

Animal Care Systems racks ship with two 6-foot-long, 3-inch hoses, each with a quick-disconnect fitting on both ends. We also supply the mating connectors used with the hoses.



Auto-Water Systems

The auto-water racks' water inlet is located at the rear top corner of the rack. The ideal rack placement along an auto-water-lined wall would be with the inlet corner against the wall. The connections are plug-and-play compatible with the industry's leading auto-water systems and are compatible with auto-flush lines.



Doorways

Animal Care Systems racks will fit through 42-inch doorways fully equipped. However, the racks also can fit through 34-inch doorways with minimal adjustments (Figure 1).




Figure 1