



Position	Lux reading
Room 316	840
Middle of room	
Top position, back of cage	60
Top position, front of cage	194
Middle position, back of cage	12
Middle position, front of cage	180
Bottom position, back of cage	6
Bottom position, front of cage	180
Room 317 middle of room	599
Top position, back of cage	46
Top position, front of cage	88
Middle position, back of cage	54
Middle position, front of cage	174
Bottom position, back of cage	16
Bottom position, front of cage	128

All light levels in the OptiMICE cages fall within the recommended guidelines according to the "Guide of Care and Use of Laboratory Animals"

Diana Carlsen
Director, Vivarium Operations
PsychoGenics
765 Old Saw Mill River Road
Tarrytown, NY 10591